

# Pedicle Subtraction Osteotomy (PSO)

## Instrument Set



### PSO Instrument Set

- Pedicle Subtraction Osteotomy (PSO) procedures provide surgeons with techniques to correct sagittal deformity
- Compatible with Momentum® Posterior Spinal Fixation System
- Curettes in both #1 and #3 cup sizes of straight, 45° up, 45° down, and 45° down with 25mm extension
- Posterior Vertebral Body Punches available in 25 & 34mm lengths
- Push button articulation, Vertebral Body Retractors in 3 different lengths
- Bayonet Nerve Root Retractors in 10 & 12mm wide tips
- Right angle L-Shaped Vertebral Body Osteotomes in 6 x 8 and 8 x 6mm
- Straight Osteotomes in 6, 10, 13mm widths

## Soft Tissue Retraction and Protection

### Adjustable Vertebral Body Retractors

- Push button articulation
- Suction holder built-in
- 3 different lengths



### Bayonet Nerve Root Retractors

- Small lip on tip designed to inhibit movement of retracted soft tissue
- Bayonet handle for visualization
- 10 & 12mm wide tips



## Vertebral Body Resection

### L-Shaped Vertebral Body Osteotomes

- Right angle osteotomes in 6 x 8 and 8 x 6mm for a simultaneous cut parallel to the nerve root and dura



### Straight Osteotomes

- Available in 6, 10, and 13mm widths to meet clinical needs



## Cancellous Bone Removal

### Curettes

- Available in both #1 and #3 cup sizes of straight, straight 45° up, 45° down, and 45° down with 25mm extension



## Fracturing the Posterior Wall

### Posterior Vertebral Body Punches

- Hook designed to grip bony posterior edge to prevent instrument movement
- Available in 25 and 34mm lengths to match patient anatomy

