

G21 V-Steady

High Viscosity Cement System



V-STEADY and
Pico Mix™

V-STEADY ADVANTAGES:

- Closed Cement Mixing System
- Long working time ~9 minutes
- 20cc Closed Cement Mixing System
- 45 sec mix time
- Waiting period after mix ~3-4 minutes
- Gun Cartridge holds up to 15cc
- Gear design for controlled cement injection
- Gun provides 1cc of cement injection per full turn
- Gun tubing adds versatility for surgeon use
- Cement consistency is sensitive to room temperature
- Check cement viscosity for desired consistency prior to injection



V-HP Gun™

Ulrich Medical USA, 18221 Edison Ave., Chesterfield, MO 63305
Phone: 866.519.0268 Fax: 636.519.0271
www.ulrichmedicalusa.com ifu.ulrichmedical.com

Over a Century
of Innovation